FlexoStep[®]

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
©Jakobs GmbH
 Version 1.0/2010

5.4.2. Inclinations:

5.4. Surface profiles and inclinations of the footboards:

Both footboards are identical.

The arrow placed on the middlebridge of the basis should always show into the direction of the exercising person.

Each footboard has two different surface structures on each side. Due to this feature the soles of feet alternatively have a specific sensory and tactile surface perception.

On one side of the footboard the surface profile is composed of cross grooves and pimples, the other side of zigzag and circular arches.

Footboards and basis are shaped asymmetrically.

This do result in 8 different inclinations in the range of 10° (slight inclination) up to 35° (strong inclination).

Generally all exercises should be started with the slightest and easiest inclination (see next page figure zigzag front down).

With increasing enhancement, with increasing safety and stability the inclination can be changed and specifically increased. This should always be adapted to the individual possibilities.

The arrow placed on the middlebridge of the basis should always show into the direction of the exercising person!

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5.3. Objectives:
Enhancement of the neuro-muscular coordination
(balance training, sensory motor skill training, equilibrium training)

Starting position stand upright

Stand upright on the FlexoStep (weight equally shared to both feet).

The device is placed on an even and skid-proof ground. The soles of your feet lie flat on the stepper.

The features of the footboard's surface (cross-grooves, pimples, zigzag, circular arches), play a minor role in this exercise and can be chosen arbitrarily.

However, you should constantly change them to stimulate the receptors in an ideal way.

Place your feet in a way that the outer knuckles are positioned about 2 to 3 cm from the FlexoStep's rotation center (both-sided notches on the footboards' surfaces).

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How to exercise

Balance both footboards into a middle position. Try to keep balance without having the footboards hit the basis at the back or front. Make sure that your heels always keep contact to the footboards!

Objective of this exercise is to keep balance. At first, secure your posture by cleaving to the ascend support or to your partner. Then loosen one and, if possible, both hands. Avoid extensive movements of the limbs or torso. Try to hold your balance by little adjustments of your legs (ankle and knee joints).

Progressions:

- Standing on one leg
- Eyes closed
- Movements of the arms, the torso and/or the head
- Additional use of the arms
- Training for specific disciplines (e. g. downhill-crouch in skiing)

b) Bending the knee

To specifically train the short parts of the calves, you should take the final position as described above.

Bend the knee-ankle forwards slightly and slowly while the heel keeps contact to the footboard. Thus you will ideally train the oftenly uncared-for muscles of calves and feet. Keep this posture for approx. 30 – 60 seconds 6 to 10 times. Relax your muscles inbetween training units by dynamic movements.

5.2.2.

Starting position sitting upright

Sitting upright (see 5.1.2.).

Place the FlexoStep so that the arrow points to the user and adjust the footboards with their small sides (pimples or zigzag) into the same direction.

This configuration at first assures the safety of the user by willfully limiting the possible extent of movement.

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The extent of movement should be raised step by step until reaching the maximum of extent in the upper ankle. If the footboards hit the back of the basis, the extent of movement is being raised due to the alteration (gyration) of the footboards and/or the FlexoStep.

However, the impact should always be soft and under control (central nervous system) to prevent injuries and overstraining.

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Progression:

- Contracture-prophylaxis
- Cartilage-gymnastics (nutrition of the ankle-cartilage)
- Enhancement of the extent of movement
- Muscle-stretching

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However, the progression should be done slowly and continually. Periodic exercising is more effective than infrequent intense training. Please prefer a long-term, progressive (“extensive”) training over short and intensive exercises.

Variations:

Keeping the final position:

At first, adjust one footboard into a save forward-position (inset to the FlexoStep).

Now put the other foot-ankle into the final position in the opposite direction (maximum pain-free extent of movement) and keep this position for about 30 – 60 seconds.

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5.2. Objective: Enhancement of flexibility (Mobilisation)

5.2.1.

Starting position stand upright

Stand upright on the FlexoStep (see 5.1.1.). Place the FlexoStep so that the arrow points to the user and adjust the footboards with their small sides (pimples or zigzag) into the same direction. This configuration at first assures the safety of the user by willfully limiting the possible extent of movement.

For progression either the footboards or the stepper can be turned around (see below: progressions / variants).

How to exercise

Set both footboards in an alternating, seesawing motion. This motion should be constant and harmonic (see 5.1.).

The extent of movement should be raised step by step until reaching the maximum of extent in the upper ankle. If the footboards hit the back of the basis, the extent of Movement is being raised due to the alteration (gyration) of the footboards and/or the FlexoStep. However, the impact should always be soft and under control (central nervous system) to prevent injuries and overstretching.

Make sure that heels and toes always keep contact to the footboards!

Progression:

- Contracture-prophylaxis
- Cartilage-gymnastics
- (nutrition of the ankle-cartilage)

- Enhancement of the extent of movement
- Endurance-training
- Muscle-stretching

The crucial parameter is the duration of the training, which depends on condition and objectives.

However, the progression should be done slowly and continually. Periodic exercising is more effective than infrequent intense training. Please prefer a long-term, progressive (“extensive“) training over short and intensive exercises.

Variations:

a) Keeping the final position:

At first, adjust one footboard into a safe forward-position (inset to the FlexoStep). Now put the other foot-ankle into the final position in the opposite direction (maximum pain-free extent of movement) and keep this position for about 30 – 60 seconds.

The slight traction that you have felt in the calve at the beginning slowly fades. Now change sides and repeat the exercise. Always train both sides.

Repeat this 6 to 10 times per training-unit. You should perform up to three units per day. Thus you stress the extension-ability of the long multi-joint muscle-parts of the calves in an ideal way. Inbetween training-units, the above described exercise leverages the blood circulation (duration: approx. 30 – 60 seconds).

This panders the physiologically healthy posture, which guarantees a perfectly harmonic movement on a natural scale combined with as much safety as possible.

While at the beginning it is better for safety and acclimatisation to look down to the feet and the FlexoStep, the view should be pointed straight-forward later (see also: Progressions). The rest of the body is relaxed according to the starting position, the respiration is aedquate.

How to exercise

Set both footboards in an alternating, seesawing motion. This motion should be constant and harmonic. The extent of movement can be raised step by step. The footboards do not hit the base on the front- nor on the backside. Thereby the movement is under permanent control of the central nervous system, which prevents injuries and overstraining. The motion goes from head to knee- and hip joints automatically ("reactively"). In the ideal case, even the spine and the shoulder girdle are being involved through the pelvic bones.

Make sure that heels and toes always keep contact to the footboards!

Progression:

- Prophylaxis of thrombosis
- Gymnastics fort he vains (enhancement of blood circulation)
- Prevention of postural hypotension
- Endurance training

The crucial parameter is the duration of the training, which depends on condition and objectives.

However, the progression should be done slowly and continually. Periodic exercising is more effective than infrequent intense training. Please prefer a long-term, progressive ("extensive") training over short and intensive exercises. Short breaks (up to 60 seconds) can raise the duration of training (interval trainig).

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To additionally enhance the body perception, perform the same exercises with asymmetrically adjusted footboards: This means that while you adjusted the crossgrooves, pimples, zigzag and circular arches equally upwards or to the front during the first exercise, you will now switch them arbitrarily. Thus you will have to develop a feeling for the extent of movement with every step, which trains your deep sensibility.

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Handling should be safe but mellow and in no way cramped, for not affecting the motion-sequence.

While at the beginning it is better for safety and acclimatisation to look down to the feet and the FlexoStep, the view should be pointed straight-forward later (see also: Progressions).

The rest of the body is relaxed according to the starting position, the respiration is adequate.

How to exercise

Set both footboards in an alternating, seesawing motion. This motion should be constant and harmonic. The extent of movement can be raised step by step.

The footboards do not hit the base on the front nor on the backside. Thereby the movement is under permanent control of the central nervous system, which prevents injuries and overstraining. The motion goes from head to knee- and hip joints automatically (“reactively”).

In the ideal case, even the spine and the shoulder girdle are being involved through the pelvic bones. Make sure that heels and toes always keep contact to the footboards!

Progression:

- Prophylaxis of thrombosis
- Gymnastics for the veins (enhancement of blood circulation)
- Prevention of postural hypotension
- Endurance training

The crucial parameter is the duration of the training, which depends on condition and objectives. However, the progression should be done slowly and continually. Periodic exercising is more effective than infrequent intense training.

Please prefer a long-term, progressive (“extensive“) training over short and intensive exercises. Short breaks (up to 60 seconds) can raise the duration of training (interval training).

While the prophylaxis of thrombosis is being performed for a short time on an hourly basis (20 – 40 seconds), exercises for the enhancement of blood circulation and stimulation of the cardiovascular system should last several minutes several times per day.

The endurance training clearly longer (up to 60 minutes) and should be scheduled daily or at least three times per week.

Variations:

To additionally enhance the body perception, perform the same exercises with asymmetrically adjusted footboards:

This means that while you adjusted the cross-grooves, pimples, zigzag and circular arches equally upwards or to the front during the first exercise, you will now switch them arbitrarily.

Thus you will have to develop a feeling for the extent of movement with every step, which trains your deep sensibility.

5.1.2.

Starting position sitting straight

Straight sitting posture on a chair or stool. The legs' weight is shared equally onto both feet. The device stands fixed on an even and skid-proof ground. The soles of your feet lie flat on the stepper.

The features of the footboard's surface (cross-grooves, pimples, zigzag, circular arches), play a minor role in this exercise and can be chosen arbitrarily. Place your feet in a way that the outer knuckles are positioned about 2 to 3 cm from the FlexoStep's rotation center (both-sided notches on the footboards' surfaces).

3. Indications

- Prophylaxis of thrombosis
- Enhancement of flexibility of ankle joints (the upper ankle joint)
- Stretching of the back muscles of both lower legs (calf muscles)
- Enhancement of endurance (prevention of cardiovascular diseases)
- Strengthening the musculature of the lower limbs: thigh and shank as well as long and short foot-musculature (in combination with weight-sleeves or elastic straps)
- Training of balance (enhancement of the neuro-muscular coordination)
Preparation for gait training

4. Warning! - security advice and protective measures

- Not suited for children under 6 years!
- Use restricted under the supervision of adults
- Weight-loading up to 120 kg.
- Keep appropriate distance to objects and peoples
- Wear light clothing for maximum range of movement
- Only exercise practice bare feet or with appropriate skid-proof shoes
- Slip hazard! (only exercise practice on skid-proof ground)

- If there is danger of falling due to dizziness, use the ascent support or grip to something
- Clarify the legs' resilience after injuries or surgery
- In case of pain or increasing existing discomforts, please seek advice from your doctor or therapist
- In case of known instabilities of the ankles or the central nervous system, please seek advice from your doctor or therapist

5. Suggestions for training

5.1 Exercises for enhancement of the blood circulation

5.1.1

Starting position stand upright

Stand upright on the FlexoStep (share your weight equally onto both feet). The device stands fixed on skid-proof-ground.

The soles of your feet are lying flat on the FlexoSteps' footboards.

The features of the footboard's surface (crossgrooves, pimples, zigzag, circular arches), play a minor role in this exercise and can be chosen arbitrarily.

Place your feet in a way that the outer knuckles are positioned about 2 to 3 cm from the FlexoStep's rotation center (both-sided notches on the footboards' surfaces). This panders the physiologically healthy posture, which guarantees a perfectly harmonic movement on a natural scale combined with as much safety as possible.

- Individuals who are exceptionally insecure as well as beginners should use the optional ascent support to lower their risk of falling and to increase their own safety. -

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Additional application of weight-sleeves or elastic straps progressively raises the training-efficiency and systematically strengthens the musculature. Defective positions as well as amyosthenia are being obviated effectively (prophylaxis).

Although the sense of balance is being stressed as well, one is able to work out one leg solely by using the stepper onesided.

The "classic" bothsided use trains the neuromuscular coordination of all areas involved. By stressing the whole bodies' balance as well as due to the posture, the overall patterns of movement of erecting areas of musculature and limbs are being trained .

The surface-sensibility of the whole-body-activity is being specially trained. The special partition and surface-structure of the footboards provides well-balanced and diversified stimulation.

The calves are being stretched gently but effectively by combining halted positions. Depending on the position of the knee joint, different areas of the calves can be trained specifically.

Even athletes or sportswoman and sportsmen will benefit from FlexoStep by using it for specific balance-training combined with typical motion-sequences.

In summary, the most important training-objectives by using FlexoStep are to improve flexibility, strenght and endurance (locally as well as overall) as well as to prevent thrombosis after injuries and surgery.

2. Target groups

In principle, FlexoStep is appropriate for people of all ages. The suggested minimum age is approx. 6 years. There is no upper age limit. Standing straight is not obligatory. The starting position depends on individual objectives. For example, FlexoStep is excellently suited for training while sitting down. Theoretically, many starting positions are possible and reasonable: the bridge-and spider-position or the half-kneeling position /stance (gunner's station).

Due to these manifold starting positions as well as the possibility to use FlexoStep onesidedly, full resilience of the respective leg is not obligatory for trainings with FlexoStep. Training can be done even with partial weight bearing, e. g. after injuries and / or surgery.

Thus, even older patients who suffer from muscular deficits or unbalancing (e. g. Dizziness) can start their training at an early stage.

Thereby, complications can be obviated and the healing process is fastened, which is supportive in rehabilitation.

For healthy individuals and athletes, FlexoStep is a real challenge in terms of sensory motor-skills. The range of use aside from usual gymnastics are prophylaxis and therapy as well as rehabilitation of certain diseases.

FlexoStep is excellent for use in physio- and occupational therapy as well as in sports.

1. Objectives

FlexoStep™ is an exercising device for improvement of motor skills (flexibility, energy, coordination and endurance), which provides manifold possibilities of application in prevention and therapy (rehabilitation).

Two footboards in light grey are directly connected to the anthracite base. Both elements of the FlexoStep are made of high-quality inherently stable and sturdy plastic which is easy to keep clean. The surface can be sanitized.

Due to the features of the footboard's surface (crossgrooves, pimples, zigzag, circular arches), sensory information is being conducted to the brain on an increased scale. This provides better perfusion of the sole as well as a higher training-efficiency.

Footboards and basis are shaped asymmetrically. The baseplate of the light but extremely stable substructure (basis) is rectangularly shaped. Thus it offers stability and ideal security to the user.

Particularly insecure or instable persons as well as beginners will find perfect additional equipment by the ascent support "Easy Step", which lowers the risk of downfall and provides a permanent secure hold while stepping on and off as well as during the whole period of exercise.

The basis-buildup features two different inclinations. Depending on the footboards' positions, the extent of movement is being altered. Also, the footboards are formed trapezoidally to adjust the extent of movement of the ankles individually to fit the user's potential as well as the training-objectives.

Constant and perpetual movement not only strenghtens the respective musculature but also enhances the blood circulation in legs and feet by altering stress and distress (thrombosis-prophylaxis).

Vulnerable individuals (older, bedridden, immobile persons / just after incidents and/or surgery) do highly benefit from the additional training with FlexoStep. This is particularly so, since the exercises can be performed sitting down, which makes it possible to treat the affected limb with caution (partial weight bearing).

The seesaw-movement makes it possible even for underperforming individuals to use FlexoStep efficiently and at an early stage; in fact, straining up- and down-movements are being avoided due to this unique principle.

Endurance grows with increasing training-time (endurance training).

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

FlexoStep®

Instruction

Jürgen Förster, BPT (B), MA
Principal of Physiotherapy
Institution of Physiotherapy
University Hospital Aachen (UKA)
Medical Faculty Technical University (RWTH)

1. Objectives (page 14)
2. Target groups (page 15)
3. Indicationen (page 16)
4. Security advice and protective measures (page 16)
5. Suggestions for training (page 16-21)
 - 5.1. Exercises for enhancement of the blood circulation (page 16-18)
 - 5.2. Objective: Enhancement of flexibility (Mobilisation) (page 19-21)
 - 5.3. Objectives: Enhancement of the neuro-muscular coordination (page 21)
 - 5.4. Surface-profiles and inclinations (page 22-23)

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